


## Sterimed Oxygen Mask

SMD 704

- Oxygen Mask
- Máscara De Oxígeno
- Máscara de oxigênio
- Masque d'oxygène



- Used for efficient oxygen administration to patient
- Made from soft, non-toxic, medical grade PVC
- Elastic strap and aluminum nose clip for proper adjustment of mask
- Multi-channel tube ensures regular supply of oxygen even if the tube is accidentally kinked
- Proximal end of tube provided with universal connector for easy connection with oxygen source

### Packing

- Non- Sterile, individually packed in a ribbon.
- Box of 20; Master Carton of 200.

Gross: 18.4kg (Adult) 740 x 425 x 425mm

SIZES	TUBE LENGTH
Adult	210cm (7ft)
Child	210cm (7ft)

## Sterimed Nasal Oxygen Cannula

Twin Bore Oxygen Delivery Set

SMD 703 A

- Nasal Oxygen Cannula Twin Bore Oxygen Delivery Set
- Canula Nasal De Oxigeno Entrega Con Doble Juego
- Nasal cânula de oxigênio Bore Twin Set entrega de oxigênio
- Canule nasale d'oxygène Bore Twin Set oxygène de livraison



SIZES	TUBE LENGTH
Adult	210cm (7ft)
Child	210cm (7ft)
Infant	210cm (7ft)

Total Length 250cm  
\* (Length can be customized)

### Packing

- Non-Sterile, individually packed in a polybag.
- Box of 25, Master Carton of 250.

Gross: 16.6kg 790 x 480 x 280mm