



Innovating For Tomorrow

# Cervi-Cath™

Cervical Ripening Balloon  
Catheter

SMD 811




Quality Since 1955

“**Sterimed Group** – your trusted healthcare partner.  
We lead with innovation and quality in providing advanced  
medical solutions. **Join us in creating a healthier world.** ”

**WHO GMP**  
Certified Company

**ISO**  
9001:2008  
13485:2012  
Certified



Visit Our Website



### Efficient Detection of Uterine and Vaginal Balloon:

Labeling of U & V on the funnel ensures precise inflation and placement during the catheter insertion.

### Advanced Flexibility with Latex Material:

The latex body of the catheter minimizes acRom risk and facilitates easy manipulation during insertion.

### Integrated Tip:

Expertly designed to prevent extra length, eliminating the risk of amniotic sac rupture.

### Double Balloon Design:

Ensures equal and steady pressure on both the internal and external cervix, minimizing displacement risk.

### Non-Pharmaceutical Option:

As a mechanical dilator, CerviCath reduces potential side effects and offers a cost-effective alternative to pharmaceutical interventions.



## Packing & Size

### Packing

- One Piece individually peelable pouch packed
- Inner Box: 1 pcs
- Second Inner Box: 10 pcs
- Master Carton: 200 pcs



Gross: 9.50kg



540 x 430 x 350mm

Size: **18 fr**

Download Our Catalogue



**STERIMED GROUP**

Regd. Office: 501, Ring Road Mall, Manglam Place, Sector-3, Rohini, New Delhi-110085 (INDIA)

Tel: +91-11-42466196 / 396 (100 Lines) | Email: info@sterimedgroup.com | Website: www.sterimedgroup.com