



Innovating For Tomorrow

Haemolock™

Postpartum Balloon Catheter

SMD 812



Quality Since 1955

“**Sterimed Group** – your trusted healthcare partner.
We lead with innovation and quality in providing advanced
medical solutions. **Join us in creating a healthier world.**”

WHO GMP
Certified Company

ISO
9001:2008
13485:2012
Certified



Visit Our Website

Haemolock - Postpartum Balloon Catheter

Open Tip:

The perforated tip allows/facilitates the drainage of the blood, clots and fluids from the uterine cavity during the tamponade process.

Rapid Instillation Method:

Ability for quick and easy insertion providing a timely response in emergency situations.



Inbuilt High Elastic Balloon:

Eliminate tiresome conventional methods and also gives maximum coverage to the uterus wall.

Non-Pharmaceutical Option:

Alternative cost effective method to address postpartum haemorrhage, minimizing the need for drugs and surgery [Allowing a less traumatic experience].

Packing & Size

Packing

- One Piece individually peelable pouch packed
- Inner Box: 1 pcs
- Master Carton: 20 pcs

Size: **24 fr**



Gross: 7.50kg



540 x 430 x 350mm

Download Our Catalogue



STERIMED GROUP

Regd. Office: 501, Ring Road Mall, Manglam Place, Sector-3, Rohini, New Delhi-110085 (INDIA)

Tel: +91-11-42466196 / 396 (100 Lines) | Email: info@sterimedgroup.com | Website: www.sterimedgroup.com