

Sterimed Microlaryngeal Tube


Oral/Nasal
SMD 765

- Tabung Microlaryngeal
- Tubo Microlaríngeo
- Tubo Microlaríngeo
- Tube Microlaryngé

Conforming to ISO 5361



PHthalate FREE



- The microlaryngeal tube's small diameter offers surgeons a clearer view and more space in the surgical field, enhancing precision during procedures.
- Ideal for patients with narrowed airways, the tube ensures optimal airflow.
- Available in DEHP-free PVC and clear PVC variants, catering to diverse patient needs while prioritizing safety.
- Used to maintain a patent airway and mechanically ventilate patients
- Made from non-toxic, non-irritant medical grade clear PVC
- Kink resistant, thin-walled, thermo-sensitive tube, softens at body temperature
- Smooth inside/outside surface for ease in intubation, extubation and suctioning
- Smooth polished Murphy eye, distal end and tapered balloon end for atraumatic intubation and extubation
- Kink resistant inflation tube to ensure patient safety during cuff inflation and deflation
- Radio opaque line and markings provided to facilitate identification of tube position
- Universal 15 mm connector at proximal end
- E.T.O Sterile and Pyrogen free. For single patient use
- Also Available in Ivory

Packing

- Sterile, Individually packed in peelable pouch / blister packed.
- Box of 10; Master Carton of 180



Gross: 10.26kg



660 x 380 x 500mm (4.0 to 6.0mm)

SIZES - CUFFED

4.0 mm, 5.0 mm, 6.0 mm