





## Sterimed Intubating Stylet

SS 753

-  intubasi stylet
-  intubación estilete
-  intubating Stylet
-  stylet d'intubation

- Used along with Endotracheal Tube to aid during intubation
- Smooth plastic coated with flexible malleable aluminum design, allows the stylet to be easily adapted to any desired shape
- Soft distal tip to reduce the possibility of intubation damage
- Smooth surface for ease in introduction and removal from ET tube
- Clean packed



### Packing

- Individually pouch packed.
- Inner Box of 10/20 pcs.
- Master Carton of 400 pcs.



Gross: 14.96kg approx.



492 x 460 x 428mm  
467 x 430 x 343mm

#### SIZES

6 Fr/Ch (2.0 mm)





10 Fr/Ch (3.3 mm)


14 Fr/Ch (4.7 mm)

## GuideX<sup>(R)</sup>

### Sterimed Introducer Bougie, Vented

SS 755

-  Pengenal Bougie, Vented
-  Introducutor Bougie, Ventilado
-  Introdutor Bougie, Ventilado
-  Introducteur Bougie, Ventilé



- Provides rigidity to endotracheal tubes to facilitate endotracheal intubation
- Ultra smooth surface for low friction
- Atraumatic angled tip for ease of intubation and to minimise risk of patient trauma
- Marking provided along the length to facilitate identification of tube position
- Angled tip for ease of intubation
- Vented lumen for oxygen feed

### Packing

- Sterile, individually packaged in easy to open sachets
- Box of 10 pcs; Master Carton of 200 Pcs.



Gross: 9.26kg approx.



790 x 345 x 350mm



Atraumatic Smooth Tip



Removable Soft Connector

#### SIZES

	6	10	15
FR/CH	6	10	15
Length	535mm	700mm	700mm