

VACPRO[®]

Suction Catheter Closed System (24 hrs / 72 hrs)

With bronchoscopy port


SMD 700 CS

- Suction Catheter Closed System (24 jam / 72 jam)
- El catéter de succión de sistema cerrado (24 hrs / 72 hrs)
- Sucção Catheter System Closed (24 horas / 72 horas)
- Des systèmes d'aspiration fermé cathéter (24 heures / 72 heures)

Conforming to ISO 8836



PHTHALATE FREE



- Used to maintain ventilation and oxygen therapy while suctioning in mechanically ventilated patients
- Minimises the risk of Ventilator Associated Pneumonia (VAP)
- Made from non-toxic, medical grade PVC
- Catheter with lateral eyes, smooth tip and depth markings for maximum patient comfort
- MDI port provides an effective, quick and convenient way of drug delivery
- Irrigation port with one way valve allows lavage and irrigation of the suction catheter
- On/off valve allows suction catheter tip to be isolated and creates an enclosed cleaning chamber. No risk of fluid going into the respiratory track
- Transparent connectors and durable sleeve facilitates catheter manipulation and prevents infiltration of blood and sputum
- 360° swivel elbow enables rotation for enhanced patient comfort and reduces torque
- Wedge for easy removal of ET Tube from the catheter.
- Flexi-mount connector for higher convenience.

SIZES						
Size	6 Fr.	8 Fr.	10 Fr.	12 Fr.	14 Fr.	16 Fr.
O.D. (mm)	2.0	2.7	3.3	4.0	4.7	5.3
I.D. (mm)	1.4	1.8	2.4	2.8	3.2	3.8
Length (mm)	300	300	600	600	600	600
Color Code	Sea Green	Blue	Black	White	Green	Orange

Packing

- One Piece individually peelable pouch packed.
- Master Carton: 50 pcs



Gross: 10.26kg



685 x 450 x 340mm