

VAPSAFE (R)

Endotracheal Tube - Subglottic Suction Line

SMD 722


- Endotracheal tube - subglottic Suction Baris
- Tubo endotraqueal - subglótica Línea de Succión
- Endotraqueal Tube - subglótica da linha de sucção
- Sonde endotrachéale - Ligne d'aspiration sous-glottique

Conforming to ISO 5361



PHthalate FREE

CE
0123



- Used to maintain a patent airway in long term mechanically ventilated patients
- Made from a special non-toxic medical grade PVC
- **Special suction port for removal of subglottic secretions using syringe/suction minimize the risk of VAP(Ventilator Associated Pneumonia)**
- The Vapsafe tube comes with a tracheal cuff
- The tracheal balloon ensures proper mucus pool formation for effective suctioning
- Extra smooth, non-adherent inside/outside surface for ease in intubation, extubation and suctioning
- Kink resistant, thin-walled, thermo-sensitive tube. Softens at body temperature
- Smooth polished murphy eye and distal end for atraumatic intubation and extubation
- Radio opaque line and markings provided to facilitate identification of tube position
- Universal 15 mm connector at proximal end
- E.T.O Sterile and Pyrogen free. For single patient use

Packing

- Sterile, Individually packed in peelable pouch / blister packed.
- Box of 25; Master Carton of 300 pcs.



Gross: 15.56kg



660 x 380 x 500mm

SIZES - CUFFED

5.0 mm, 5.5 mm, 6.0 mm, 6.5 mm, 7.0 mm, 7.5 mm, 8.0 mm, 8.5 mm