



## Sterimed Yankaur Suction Set

With Standard Handle | SMD 712

With Crown Tip Handle | SMD 709

- Yankaur Suction Set With Standard Handle
- Juego De Aspiración Yankaur Con Asa Standard
- Yankaur Conjunto de Sucção Com Punho Padrão
- Yankaur aspiration avec poignée standard

- Used for removal of secretion and blood pass during surgical procedures
- Made from non-toxic, non-irritant PVC
- Flexible kink resistant tube provides uninterrupted flow during suction
- Soft connectors at both ends of the tube for easy connection
- Available in Standard / Round Tip Handle
- Suitable for use with or without vent control facility

### Packing

- Sterile, individually packed in peelable pouch pack.
- Box of 10; Master Carton of 100



Gross: 17.5kg



740 x 425 x 425mm



PHTHALATE FREE



SIZES	
Tube Length	2500mm (2.5mtr.)
Tube Diameter	9mm



## Stericiser

Breathing Exerciser

SS 1031

- Berolahraga Pernapasan
- Ejercitador Respiratorio
- respiração Exerciser
- respiration exerciseur

- Used as a breathing exerciser to increase lung capacity in patients with decreased lung function
- Made from rigid, non-toxic medical grade PVC
- Specially designed with 3 stage chambers to provide more efficient exercise to patient at a step up basis
- Device comprises a base and a central part with a 3 stage chamber, each containing spheres of different sizes, weight and colour
- Provides lung exercise on a step-up basis, as the spheres rise proportionately to the breathing effort exerted by the patient
- Can be dis-assembled into parts for cleaning and disinfection
- Non-sterile

### Packing

- Sterile / Disposable / Individual Printed Box Packing.
- Box of 1/Master Box of 50.



Gross: 10.5kg approx.



660 x 380 x 420mm